


Seurat scRNA pipeline

Mingbo Cheng


13/12/2020

Idea

Data Producing





Visualization



How it works

Data Producing



meta.data

clustering results
annotation
.....

reductions


PCA
UMAP
.....

tools

DE analysis
GO analysis
Pathways analysis
.....




Visualization




How it works


Data Producing




How it works




How it works



How it works




How it works



How it works

Data Producing



conf/config.R


data_factory.R

external marks etc.

Basic info: Species, Organ etc.


```
### -----Tinitail info-----  
PROJECT = "Mouse heart Gli&CD45 project" ## set project name  
ORGAN = 'Heart' #For external annotation. Options: Blood,  
SPECIES = "Mouse" #For external annotation. Options: Human,  
MCA_NAME = "Neonatal-Heart" #For MCA annotation. Options: check  
  
# filtering params when create seurat object  
MINCELLS = 5  
MINGENES = 50  
  
### ----- Data SRC-----  
ANNOTATION_EXTERNAL_FILE = "external/Human_and_mouse_cell_markers-M  
  
data_src = c(  
  NK1_Gli1_IRI = "data/summed_mtx/sum_NK1_Gli1_IRI/",  
  NK2_CD45_IRI = "data/summed_mtx/sum_NK2_CD45_IRI/",  
  NK3_Gli1_Shram = "data/summed_mtx/sum_NK3_Gli1_Shram/",  
  NK4_CD45_Shram = "data/summed_mtx/sum_NK4_CD45_Shram/"  
)  
  
##----- SET REPLICATE GROUP -----  
stage_lst = c(  
  NK1_Gli1_IRI = "IRI",  
  NK2_CD45_IRI = "IRI",  
  NK3_Gli1_Shram = "Shram",  
  NK4_CD45_Shram = "Shram"
```

How it works



How it works

Data Producing



conf/config.R



data_factory.R




external marks etc.

Basic info: Species, Organ etc.


Input location: 10x/h5

```
### -----Initial Info-----  
PROJECT = "Mouse heart Gli&CD45 project" ## set project name  
ORGAN = 'Heart' #For external annotation. Options: Blood, Heart  
SPECIES = "Mouse" #For external annotation. Options: Human, Mouse  
MCA_NAME = "Neonatal-Heart" #For MCA annotation. Options: check http  
  
# filtering params when create seurat object  
MINCELLS = 5  
MINGENES = 50  
  
### ----- Data SRC-----  
ANNOTATION_EXTERNAL_FILE = "external/Human_and_mouse_cell_markers-Markers  
  
data_src = c(  
  NK1_Gli1_IRI = "data/summed_mtx/sum_NK1_Gli1_IRI/",  
  NK2_CD45_IRI = "data/summed_mtx/sum_NK2_CD45_IRI/",  
  NK3_Gli1_Shram = "data/summed_mtx/sum_NK3_Gli1_Shram/",  
  NK4_CD45_Shram = "data/summed_mtx/sum_NK4_CD45_Shram/"  
)  
  
##----- SET REPLICATE GROUP -----  
stage_lst = c(  
  NK1_Gli1_IRI = "IRI",  
  NK2_CD45_IRI = "IRI",  
  NK3_Gli1_Shram = "Sham",  
  NK4_CD45_Shram = "Sham"  
)
```


How it works




How it works




How it works



How it works



How it works





!!!Know every step

```
## 0. omit,      1. calc & save,      2. load      3.
conf = c(
  scrna_rawdata      = 1, ## read count matrix and merge samples to a Seurat obj
  scrna_filter       = 1, ## filter nFeatureRNA and nCountRNA
  scrna_preprocess   = 1, ## Normalize & FindvariableFeatures and ScaleData
  scrna_cellcycle    = 1, ## Cell cycle scoring
  scrna_cycleRegressOut = 1, ## Regress out cell cycle effects
  scrna_CCmitoRegressOut = 1, ## Regress out mito & cell cycle
  scrna_mitoRegressOut = 0, ## Regress out mito only
  scrna_riboRegressOut = 0, ## Regress out ribo only
  scrna_CCriboRegressOut = 0, ## Regress out ribo & cell cycle
  scrna_mitoRiboRegressOut = 0, ## Regress out ribo&mito
  scrna_CCmitoRiboRegressOut = 0, ## Regress out ribo&mito & cell cycle
  scrna_integration  = 1, ## Integrate samples using Seurat 3
  scrna_ScaleIntegration = 1, ## ScaleData&PCA and UMAP
  scrna_batchclustering = 1, ## clustering with resolution from 0.1 to 0.8
  scrna_batch_markersgenes = 1, ## Marker Genes for clusters with different resolutions
  scrna_clustering     = 1, ## Set seurat_clusters or re-calculate
  scrna_clusterwise_xcell = 1, ## remove cells of each cluster according distinct criterion
  scrna_cluster_annotation = 1, ## Annotate clusters according to `cluster_annotation`
  scrna_del_mitogenes   = 0, ## !!!DANGEROUS, once deleted, never recovered!!!
  scrna_markersgenes    = 1, ## markersgenes for seurat_clusters
  scrna_genesortR       = 1, ## genesortR analysis
  scrna_go              = 1, ## Gene Ontology analysis
  scrna_kegg            = 1, ## kegg enrichment analysis
  scrna_reactome        = 1, ## reactome enrichment analysis
  scrna_hallmark        = 1, ## hallmark enrichment analysis
  scrna_fishertest_clusters = 1, ## fisher test for clusters and stages
  scrna_MCAannotate     = 1, ## scMCA annotation celltypes
  scrna_ExternalAnnotation = 1, ## Annotation from given databases(tsv)
  scrna_dego_name       = 1, ## DE & GO between samples
  scrna_dego_stage      = 1, ## DE & GO between stages
  scrna_dego_stage_vsRest = 1, ## DE & GO between one stage and all Rest
  scrna_pathway_name     = 1, ## samples comparison KEGG&Reactome&hallmark
  scrna_pathway_stage    = 1, ## stages comparison KEGG&Reactome&hallmark
  scrna_pathway_stage_vsRest = 1, ## stages vsRest comparison KEGG&Reactome&hallmark
  scrna_clusterwise_xcell = 0, ## keep cells for each cluster according to mito&ribo
  scrna_fishertest_clusters = 0, ## fisher test for clusters and stages
  scrna_merge_clusters   = 0, ## merge clusters
  scrna_remove_clusters  = 0, ## remove clusters
  scrna_remove_recluster = 0, ## remove clusters and recluster with default resolution
  scrna_markersgenes     = 0, ## markersgenes for seurat_clusters
  scrna_go               = 0, ## Gene Ontology analysis
  scrna_dego_name        = 0, ## DE & GO between samples
  scrna_dego_stage       = 0) ## GO down for mark genes
```

How it works

Data Producing



meta.data

clustering results
annotation
.....

reductions


PCA
UMAP
.....

tools

DE analysis
GO analysis
Pathways analysis
.....




Visualization




How it works

Visualization



How it works



How it works

```
#!/bin/bash
RED='\033[0;31m'
NC='\033[0m' # No Color
```

```
FUNCS=(
    QC
    DEs
    Clusters
    # DEGO
    ## EXT_MARKERS
    DEGO_1v1
    DEGO_stage
)
```

```
#!!!!!!!!!!!!-----clusters to choose-----
# In general, we choose seurat_clusters,
# If you are using removed or merged clusters,
# choose the following:
# seurat_clusters
# merged_clusters
# removed_clusters
# remove_recluster
```

```
#cluster="removed_clusters"
#cluster="remove_recluster"
#cluster="merged_clusters"
#cluster="annotation"
#cluster="singleton"
```

```
cluster="seurat_clusters"
```

```
#!!!!!!!!!!!!-----
```


```
echo -e "Use cluster slot ${RED} $cluster ${NC}"
mkdir -p report/data
```

- conf/config.R
- excels generated
- Seurat Robjects


run.sh

- 1_quality_report.Rmd
- 2_clustering.Rmd
- 2_clusters_DEs.Rmd
- 3_DE_GO-analysis.Rmd
- 3_external_markers.Rmd
- make_report.R
- singleton_clustering.Rmd
- template/
 - DE-GO-1v1.template
 - DE-GO-stagesVS.template
 - index.template

Visualization



How it works




Most Frequent questions

- Is it possible to do further analysis given an existed Seurat Object?

Most Frequent questions

- Is it possible to do further analysis given an existed Seurat Object?



Most Frequent questions

- Is it possible to do further analysis given an Seurat Object?

```
## 0. omit,      1. calc & save,      2.  
conf = c(  
  scrna_rawdata           = 0,  
  scrna_filter            = 0,  
  scrna_preprocess        = 0,  
  scrna_cellcycle         = 0,  
  scrna_cycleRegressOut   = 0,  
  scrna_CCmitoRegressOut  = 0,  
  scrna_mitoRegressOut    = 0,  
  scrna_riboRegressOut    = 0,  
  scrna_CCriboRegressOut  = 0,  
  scrna_mitoRiboRegressOut = 0,  
  scrna_CCmitoRiboRegressOut = 0,  
  scrna_integration       = 0,  
  scrna_ScaleIntegration  = 0,  
  scrna_batchclustering   = 0,  
  scrna_batch markergenes = 0,  
  scrna_clustering        = 2,  
  scrna_clusterwise_xcell = 1,  
  scrna_cluster_annotation = 1,  
  scrna_del_mitogenes     = 0,  
  scrna markergenes       = 1,
```

Most Frequent questions


- Is it possible to do further analysis given an Seurat Object?

!!! umap reduction name

```
## 0. omit,      1. calc & save,      2.
conf = c(
  scrna_rawdata          = 0,
  scrna_filter           = 0,
  scrna_preprocess       = 0,
  scrna_cellcycle        = 0,
  scrna_cycleRegressOut  = 0,
  scrna_CCmitoRegressOut = 0,
  scrna_mitoRegressOut   = 0,
  scrna_riboRegressOut   = 0,
  scrna_CCriboRegressOut = 0,
  scrna_mitoRiboRegressOut = 0,
  scrna_CCmitoRiboRegressOut = 0,
  scrna_integration      = 0,
  scrna_ScaleIntegration = 0,
  scrna_batchclustering  = 0,
  scrna_batch markergenes = 0,
  scrna_clustering       = 2,
  scrna_clusterwise_xcell = 1,
  scrna_cluster_annotation = 1,
  scrna_del_mitogenes     = 0,
  scrna markergenes       = 1,
```

Most Frequent questions

- Is it possible to do further analysis given an Seurat Object?
- Why there are so many warnings especially for GO & pathway analysis?



Most Frequent questions

- Is it possible to do further analysis given an Seurat Object?
- Why there are so many warnings especially for GO & pathway analysis?
- Failed to run visualization part: cannot find 'INTE_UMAP'

Most Frequent questions

- Is it possible to do further analysis given an Seurat Object?
- Why there are so many warnings especially for GO & pathway analysis?
- Failed to run visualization part: cannot find 'INTE_UMAP'
- Why there are so many Robjects generated?

Most Frequent questions

- Is it possible to do further analysis given an Seurat Object?
- Why there are so many warnings especially for GO & pathway analysis?
- Failed to run visualization part: cannot find 'INTE_UMA'?
- Why there are so many Robjects generated?

```
180M Oct 22 09:21 scrna_RegressOutAll.Rds
121M Nov 10 16:20 scrna_go.Rds
121M Nov 10 16:20 scrna_fishertest_clusters.Rds
121M Nov 10 16:21 scrna_MCAannotate.Rds
121M Nov 10 16:22 scrna_ExternalAnnotation.Rds
198M Nov 10 16:36 scrna_dego_name.Rds
218M Nov 10 16:40 scrna_kegg.Rds
220M Nov 10 16:42 scrna_reactome.Rds
220M Nov 10 16:43 scrna_hallmark.Rds
220M Nov 10 16:47 scrna_dego_stage.Rds
221M Nov 10 16:49 scrna_pathway_stage.Rds
224M Nov 10 16:53 scrna_pathway_name.Rds
2.6M Nov 10 16:54 all_de_list.Rds
262M Nov 10 16:59 scrna_dego_stage_vsRest.Rds
263M Nov 10 17:02 scrna_pathway_stage_vsRest.Rds
740M Nov 11 17:49 scrna_rawdata.Rds
13M Nov 11 17:50 scrna_filter.Rds
166M Nov 11 17:50 scrna_preprocess.Rds
20M Nov 11 17:50 scrna_cellcycle.Rds
175M Nov 11 17:50 scrna_cycleRegressOut.Rds
176M Nov 11 17:51 scrna_CCmitoRegressOut.Rds
177M Nov 11 17:51 scrna_mitoRegressOut.Rds
178M Nov 11 17:51 scrna_riboRegressOut.Rds
180M Nov 11 17:52 scrna_CCriboRegressOut.Rds
181M Nov 11 17:52 scrna_mitoRiboRegressOut.Rds
183M Nov 11 17:52 scrna_CCmitoRiboRegressOut.Rds
28M Nov 11 17:53 scrna_integration.Rds
45M Nov 11 17:53 scrna_ScaleIntegration.Rds
45M Nov 11 17:53 scrna_batchclustering.Rds
47M Nov 11 18:05 scrna_batch_markergenes.Rds
47M Nov 11 18:05 scrna_clustering.Rds
38M Nov 11 18:05 scrna_clusterwise_xcell.Rds
38M Nov 11 18:06 scrna_markergenes.Rds
46M Nov 30 20:37 scrna_del_mitogenes.Rds
```

TODO

- Merge code from Tiago
- Add harmony integration
- NABA geneset score
- Add KEGG/Reactome/hallmark visualization
- scHCL for human cell annotation
- Integrated with ligand receptor analysis?

Dummy example

- our hpc : 134.130.18.27
 - module add scRNA/1.0.3
 - /data/scRNA/scrna_seurat_pipeline_demo

Thanks

Q&A